

Detection of cytokine and Nrf2 signaling pathway component expression in THP-1 cells

Quantitative real-time reverse-transcription polymerase chain reaction (qRT-PCR)

AUTHORED BY:	DATE:
Cordula Hirsch	16.01.2014

REVIEWED BY:	DATE:
Harald Krug	10.04.2014

APPROVED BY:	DATE:

DOCUMENT HISTORY

Effective Date	Date Revision Required	Supersedes
15.02.2014	DD/MM/YYYY	DD/MM/YYYY

Version	Approval Date	Description of the Change		Author / Changed by
1.0	DD/MM/YYYY	All	Initial Document	Cordula Hirsch

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		1/20

Table of Content

1	Introduction.....	3
2	Principle of the Method	3
3	Applicability and Limitations	3
4	Related Documents	4
5	Equipment and Reagents	4
5.1	Equipment	4
5.2	Reagents	5
5.3	Reagent Preparation	6
5.3.1	Complete cell culture medium	6
5.3.2	PMA stock solution.....	6
5.3.3	NaOH	6
5.3.4	Pluronic F-127.....	6
5.3.5	LPS	7
6	Procedure	7
6.1	General remarks.....	7
6.2	Flow chart 1	8
6.3	Cell seeding.....	9
6.3.1	Cell culture.....	9
6.3.2	Cell seeding into 24-well plate	9
6.4	Cell treatment.....	10
6.4.1	Dilution of nanomaterials.....	10
6.4.2	Dilution of LPS (chemical positive control).....	12
6.4.3	Application of stimuli.....	12
6.5	RNA isolation	13
6.6	RNA measurement	15
6.7	cDNA synthesis using iScript™ (BioRad).....	16
6.8	Real-time PCR using iQ™ SYBR® Green Supermix and iCycler	16
6.9	Data evaluation	18
7	Quality Control, Quality Assurance, Acceptance Criteria.....	18
8	Health and Safety Warnings, Cautions and Waste Treatment.....	19
9	Abbreviations	19
10	References.....	20

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		2/20

1 Introduction

Oxidative stress reactions are likely to be induced upon nanomaterial (NM) – cell interactions (e.g. Donaldson et al. 2004; Schins and Knaapen, 2007). An imbalance between the production of reactive oxygen species (ROS) and their degradation is defined as oxidative stress (e.g. Limbach et al. 2007). Cells initiate defense mechanisms against ROS overproduction such as the Nrf2 (nuclear factor erythroid 2-related factor 2) signaling pathway to assure cell survival. The induction of Nrf2 expression as well as expression of its target genes can be monitored by qRT-PCR analysis.

In case of failure of the antioxidant defense system, ROS may lead to damage of DNA, proteins and lipids. Such injuries may subsequently lead to genotoxicity, inflammation and cell death. In general inflammation describes a systemic and complex reaction of the body to harmful stimuli, as e.g. pathogens or irritants. This process involves (among others) the production of different cytokines by different cell types to allow for a coordinated defense reaction of the body. Apart from analyzing cytokine release using the ELISA technique (enzyme-linked immunosorbent assay; see SOP “Detection of cytokine release in THP-1 cells”) the transcription level of the respective cytokines can be studied by qRT-PCR analysis.

2 Principle of the Method

Quantitative real-time reverse transcription polymerase chain reaction (qRT-PCR) is used to amplify and simultaneously quantify a DNA molecule (gene) of interest. The procedure relies on the general PCR principle which will not be further detailed here (for basic textbook knowledge see e.g. “Molekulare Genetik, Rolf Knippers”).

To assess the amount of an actively transcribed gene of interest its level of messenger RNA (mRNA) is analyzed. Therefore total RNA is isolated and purified. Using reverse transcriptase copy DNA (cDNA) is synthesized and serves as the template for PCR analysis.

Inflammatory reactions are analyzed using specific primers for human interleukin 8 (IL-8) and human tumor necrosis factor alpha (TNF- α). To assess Nrf2-signaling specific primers for human Nrf2 and human heme oxygenase 1 (HO-1) are used. Glyceraldehyde 3-phosphate dehydrogenase (GAPDH) is analyzed in parallel as the housekeeping gene. For PCR reaction and DNA detection an iCycler (BioRad) and Sybr Green are used, respectively. Gene expression levels are determined using the $\Delta\Delta CT$ method.

3 Applicability and Limitations

Cytokine expression and release is cell type dependent. Not all cell types release cytokines and not all cytokines are released by one cell type. This SOP specifically addresses the analysis of IL-8 and TNF- α expression in THP-1 cells. To be able to directly compare gene regulation on the mRNA level and protein release (assessed by ELISA, described in SOP “Detection of cytokine release in THP-1 cells”) we harvest both biomolecules from the same sample. This necessitates optimized culture conditions as described in chapter 6 “Procedure”. We use the RNeasy Micro kit (Qiagen) for RNA extraction and

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		3/20

purification, the iScript cDNA synthesis kit (BioRad) for cDNA synthesis and SYBR Green to detect DNA amplification in an iCycler (BioRad). Any other mRNA extraction method and cDNA synthesis procedure revealing comparable purity and yield of mRNA can be used instead.

4 Related Documents

Table 1: Documents needed to proceed according to this SOP and additional NM-related interference control protocols.

Document ID	Document Title
cell culture_THP-1	<i>Culturing of THP-1 cells</i>
	<i>RNeasy® Micro Handbook - Qiagen</i>
qRT-PCR_efficiency	<i>Efficiency testing for qRT-PCR - $\Delta\Delta CT$ method for qRT-PCR data analysis</i>
M_NM suspension_metal oxides	<i>Suspending and diluting Nanomaterials – Metal oxides and NM purchased as monodisperse suspensions</i>
M_NM suspension_carbon based	<i>Suspending and diluting Nanomaterials – Carbon based nanomaterials</i>

5 Equipment and Reagents

5.1 Equipment

- Centrifuge (for cell pelleting; able to run 15 ml as well as 50 ml tubes at 200 x g)
- Conical tubes (15 ml and 50 ml; polypropylene or polystyrene; e.g. from Falcon)
- Flat bottom 24-well cell culture plates
- Hemocytometer
- iCycler (or comparable real-time PCR device)
- Laminar flow cabinet (biological hazard standard)
- Light microscope (for cell counting and cell observation)
- Microcentrifuge (with rotor for 2 ml tubes)
- Micro-pipettes (0.5 μ l – 1 ml)
- Microreaction tubes (1.5 ml; e.g. from Eppendorf)
- optically clear heat seal (e.g. BioRad)
- qRT-PCR suitable 96-well plates (e.g. Multiplate™ low-profile 96-well unskirted PCR plates, BioRad)
- RNase-free pipet tips
- spectrophotometer (to assess RNA concentration and purity; e.g. Nanodrop ND-1000)
- Vortex®

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		4/20

5.2 Reagents

For cell culturing and differentiation:

- Dimethyl sulfoxide (DMSO) [CAS number: 67-68-5]
- Fetal Calf Serum (FCS)
- L-glutamine
- Lipopolysaccharides (LPS) from *Salmonella enterica* [Sigma L7770; no CAS number]
- Neomycin¹⁾
- Penicillin¹⁾
- Phorbol 12-myristate 13-acetate (PMA) [CAS number: 16561-29-8]
- **Note: Carcinogenic! Handle with special care! Special waste removal** (see chapter 8)
- Phosphate buffered saline (PBS)
- Roswell Park Memorial Institute medium (RPMI-1640)
- Streptomycin¹⁾

¹⁾ bought as a 100x concentrated mixture of Penicillin, Streptomycin and Neomycin (PSN) e.g. from Gibco.

Additionally necessary to dilute carbon based NM:

- 10x concentrated RPMI-1640
- Sodium bicarbonate solution, 7.5% (NaHCO₃) [CAS-number: 144-55-8]

For RNA isolation, cDNA synthesis and qRT-PCR:

- DNase I [Qiagen #79254]
- Ethanol (70% and 100%) [CAS number: 64-17-5]
- iQTM SYBR[®] Green Supermix [BioRad#170-8882]
- iScriptTM cDNA Synthesis Kit [BioRad #170-8891]
- Nuclease-free ddH₂O
- QIAshredder homogenizer [Qiagen #79656]
- RNeasy[®] Micro Kit [Qiagen #74004]
- β-mercaptoethanol [CAS number: 60-24-2]
- specific primer pairs:

name	gene	sequence	amplicon length	annealing temp.
hIL8_sense hIL8_antisense	human interleukin 8	5'-tct gca gct ctg tgt gaa gg-3' 5'-aat ttc tgt gtt ggc gca gt-3'	153 bp	60°C
hTNFα_sense hTNFα_antisense	human tumor necrosis factor alpha	5'-ccg tct cct acc aga cca ag-3' 5'-ctg agt cgg tca ccc ttc tc-3'	148 bp	60°C
hNrf2_sense hNrf2_antisense	human nuclear factor erythroid 2-related factor2	5'- aac cac cct gaa agc aca gc-3' 5'- tga aat gcc gga gtc aga atc-3'	51 bp	60°C
hHO-1_sense hHO-1_antisense	human heme oxygenase-1	5'- atg gcc tcc ctg tac cac atc-3' 5'- tgt tgc gct caa tct cct cct-3'	55 bp	60°C
GAPDH_sense GAPDH_antisense	human glyceraldehyde 3-phosphate dehydrogenase	5'-agt cag ccg cat ctt ctt tt-3' 5'-cca ata cga cca aat ccg ttg-3'	97 bp	60°C

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		5/20

For nanomaterial (NM) suspension and dilution:

- Pluronic F-127 [CAS number: 9003-11-6]

For waste treatment:

- HCl (smoking) [CAS number: 7647-01-0]
Note: Corrosive and Irritating! Handle with special care! (see chapter 8)
- NaOH [CAS number: 1310-73-2]
Note: Corrosive! Handle with special care! (see chapter 8)

5.3 Reagent Preparation

5.3.1 Complete cell culture medium

Basic medium:

- RPMI-1640

supplemented with:

- 10% FCS
- 1x PSN, which results in final concentrations of:
 - 50 µg/ml Penicillin
 - 50 µg/ml Streptomycin
 - 100 µg/ml Neomycin
- 0.2 mg/ml L-glutamine

5.3.2 PMA stock solution

Prepare a 1 mM stock of PMA in DMSO. Therefore resuspend the 1 mg (standard packaging size) PMA powder in 1.62 ml DMSO. Aliquote and freeze at -20°C. Can be stored for years.

Note: Carcinogenic! Handle with special care! Special waste removal (see chapter 8)

5.3.3 NaOH

Prepare a 5 M solution NaOH for PMA waste treatment.

- Dissolve 200 g NaOH pellets in 1 l ddH₂O.

Note: Be careful, exothermic reaction, gets HOT. NaOH is corrosive, wear protective clothing (especially eye protection).

5.3.4 Pluronic F-127

Stock:

- 160 ppm in ddH₂O: 160 µg/ml (=16 mg/100 ml)

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		6/20

5.3.5 LPS

Stock:

- 1 mg/ml in PBS: reconstitute the whole vial (1 mg) in 1 ml of sterile PBS
- Freeze this stock in 100 µl aliquots at -80°C.
- Can be stored for years.

Single use working stock:

- Dilute the 1 mg/ml stock solution to 100 µg/ml in PBS (1:10).
- Freeze this 100 µg/ml working stock as 50 µl **single use** aliquots at -80°C.
- Never re-freeze after thawing!

6 Procedure

6.1 General remarks

Well size and cell numbers are optimized to allow mRNA and protein measurements from the same sample. Cells (lysed to obtain mRNA) and supernatants (containing proteins) are harvested after 3, 8 and 24 hours of treatment. For technical reasons a separate 24-well plate for each time point is used.

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		7/20

6.2 Flow chart 1




Figure 1: Brief outline of the workflow: from cell seeding to qRT-PCR performance.

*RLT-working buffer is specific for the RNeasy Micro kit (Qiagen). Use the sampling buffer suitable for your RNA extraction method. Samples can be stored in RLT-working buffer at -20°C for several weeks or even months.

**At this step the supernatant of the cells can be harvested for cytokine release analysis by ELISA (see SOP "Detection of cytokine release in THP-1 cells").

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		8/20

6.3 Cell seeding

6.3.1 Cell culture

THP-1 cells are grown in T75 cell culture flasks in a total volume of 20 ml of complete cell culture medium. They are kept at 37°C, 5% CO₂ in humidified air in an incubator (standard growth conditions according to SOP “Culturing and differentiating THP-1 cells”).

6.3.2 Cell seeding into 24-well plate

- Three days (72 h) prior to experimental start harvest and count cells as described in SOP “Culturing and differentiating THP-1 cells”.
- Seed 2.5×10^5 cells in 500 μ l complete cell culture medium containing 200 nM PMA per well into a 24-well cell culture plate.
- For one 24-well plate (see Figure 2) 2.5×10^6 cells are suspended in 5 ml complete PMA containing cell culture medium (5×10^5 cells/ml). To assess 3 time points, 3 plates are necessary (7.5×10^6 cells in 15 ml).

Note: PMA is diluted 1:5000 from the 1 mM stock (2 μ l/10 ml medium).

- Using a 1 ml micro-pipette 500 μ l of this cell suspension are distributed into each of the green wells (B1 to B3 and C1 to C6, Figure 2).




Figure 2: Cell seeding into a 24-well plate.

Cells are seeded at a density of 2.5×10^5 cells per well in 500 μ l complete cell culture medium (containing 200 nM PMA) into each of the green wells. Black wells receive 500 μ l complete cell culture medium each.

- Remaining wells (labeled in black in Figure 2) receive 500 μ l complete cell culture medium only.
- Differentiate cells for three days (72 hours) in a humidified incubator at standard growth conditions.

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		9/20

6.4 Cell treatment

6.4.1 Dilution of nanomaterials

For this SOP we distinguish two types of nanomaterials (NM) according to their solvent, suspension properties and highest concentrations used in the assay. See also respective related documents (3).

- (1) Metal oxide NM, Polystyrene beads and all NM delivered as monodisperse suspensions by the supplier: solvent either determined by the supplier or ddH₂O; sub-diluted in ddH₂O; highest concentration in assay 100 µg/ml
- (2) Carbon based NM: suspended and sub-diluted in 160 ppm Pluronic F-127; highest concentration in assay 80 µg/ml

Volumes given in the following dilution schemes are enough for three 24-well plates.

Note: "Mixing" in the context of diluting NMs means, the solvent containing tube is put on a continuously shaking Vortex® and the previous sub-dilution (or stock suspension, respectively) is put dropwise into the shaking solvent. The resulting suspension stays on the Vortex® for additional 3 seconds before proceeding with the next sub-dilution.

(1) Metal oxide NM:

Prepare serial sub-dilutions of the stock suspension (1 mg/ml) in ddH₂O:

- Label six microreaction tubes (1.5 ml total volume) with 1 to 6 (relates to steps 1-6 below).
 - Add 1 ml NM stock suspension to tube no. 1.
 - Add 350 µl ddH₂O to tubes no. 2, 4, 5 and 6.
 - Add 390 µl ddH₂O to tube 3.
1. 1 ml NM stock suspension in ddH₂O → 1 mg/ml (1)
 2. 350 µl of 1 mg/ml stock suspension (1) are mixed with 350 µl of ddH₂O → 500 µg/ml (2)
 3. 260 µl of 500 µg/ml (2) are mixed with 390 µl ddH₂O → 200 µg/ml (3)
 4. 350 µl of 200 µg/ml (3) are mixed with 350 µl ddH₂O → 100 µg/ml (4)
 5. 350 µl of 100 µg/ml (4) are mixed with 350 µl ddH₂O → 50 µg/ml (5)
 6. 390 µl ddH₂O → solvent control (6)

Preparation of final dilutions:

- Label six conical tubes (15 ml total volume) as follows:
 1. 100 µg/ml
 2. 50 µg/ml
 3. 20 µg/ml
 4. 10 µg/ml
 5. 5 µg/ml
 6. Solvent control
- Add 1.8 ml complete cell culture medium to each tube.
- Mix on the Vortex® with 200 µl of the respective NM sub-dilutions or the solvent (ddH₂O):

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		10/20

1. 200 µl of the stock suspension (1 mg/ml) are mixed with 1.8 ml medium
→ 100 µg/ml (1)
2. 200 µl of the 500 µg/ml sub-dilution are mixed with 1.8 ml medium → 50 µg/ml (2)
3. 200 µl of the 200 µg/ml sub-dilution are mixed with 1.8 ml medium → 20 µg/ml (3)
4. 200 µl of the 100 µg/ml sub-dilution are mixed with 1.8 ml medium → 10 µg/ml (4)
5. 200 µl of the 50 µg/ml sub-dilution are mixed with 1.8 ml medium → 5 µg/ml (5)
6. 200 µl of ddH₂O (solvent) are mixed with 1.8 ml medium → solvent control (6)

(2) Carbon based NM:

Prepare serial sub-dilutions of the stock suspension (500 µg/ml) in 160 ppm Pluronic F-127:

- Label six microreaction tubes (1.5 ml total volume) with 1 to 6 (relates to steps 1-6 below).
 - Add 1 ml NM stock suspension to tube no. 1.
 - Add 500 µl 160 ppm Pluronic F-127 to tubes 2 to 6.
1. 1 ml NM stock suspension in 160 ppm Pluronic → 500 µg/ml (1)
 2. 500 µl of 500 µg/ml stock suspension (1) are mixed with 500 µl of Pluronic F-127
→ 250 µg/ml (2)
 3. 500 µl of 250 µg/ml (2) are mixed with 500 µl Pluronic F-127 → 125 µg/ml (3)
 4. 500 µl of 125 µg/ml (3) are mixed with 500 µl Pluronic F-127 → 62.5 µg/ml (4)
 5. 500 µl of 62.5 µg/ml (4) are mixed with 500 µl Pluronic F-127 → 31.25 µg/ml (5)
 6. 500 µl 160 ppm Pluronic F-127 → solvent control (6)

Preparation of final dilutions:

- Prepare the appropriate dilution of a 10x concentrated medium stock as follows. This mixture (A) is used in all following steps for the preparation of the final NM concentrations. Mixing NM sub-dilutions with (A) will result in 1x concentrated medium containing the correct concentrations of all supplements and the respective NM concentrations.

Reagent	Volume
10x RPMI	1.7 ml
100x PSN	170 µl
100x L-Glutamine	170 µl
7.5% NaHCO ₃	450 µl
100% FCS	1.7 ml
ddH ₂ O	10 ml

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		11/20

- Label six conical tubes (15 ml total volume) as follows:
 1. 80 µg/ml
 2. 40 µg/ml
 3. 20 µg/ml
 4. 10 µg/ml
 5. 5 µg/ml
 6. Solvent control
- Add 2.1 ml (A) to each tube.
- Mix on the Vortex® with 400 µl of the respective NM sub-dilutions or the solvent (160 ppm Pluronic F-127):
 1. 400 µl of the stock suspension (500 µg/ml) are mixed with 2.1 ml medium → 80 µg/ml (1)
 2. 400 µl of the 250 µg/ml sub-dilution are mixed with 2.1 ml medium → 40 µg/ml (2)
 3. 400 µl of the 125 µg/ml sub-dilution are mixed with 2.1 ml medium → 20 µg/ml (3)
 4. 400 µl of the 62.5 µg/ml sub-dilution are mixed with 2.1 ml medium → 10 µg/ml (4)
 5. 400 µl of the 31.25 µg/ml sub-dilution are mixed with 2.1 ml medium → 5 µg/ml (5)
 - 400 µl of 160 ppm Pluronic F-127 (solvent) are mixed with 2.1 ml medium → solvent control (6)

6.4.2 Dilution of LPS (chemical positive control)

Prepare a 10 µg/ml sub-dilution of the LPS working stock (100 µg/ml) in PBS:

- mix 27 µl of PBS with 3 µl of the working stock.

Prepare the final concentrations in complete cell culture medium as follows:

- 1 µg/ml: 2 ml medium + 20 µl working stock (100 µg/ml)
- 0.1 µg/ml: 2 ml medium + 20 µl sub-dilution (10 µg/ml)

Apply NM as well as LPS as described below.

6.4.3 Application of stimuli

Note: All NM dilutions have to be vortexed directly before application to the cells.

Three days (72 hours) after cell seeding!

- Remove PMA-containing medium. **Toxic! Discard this supernatant separately! See Chapter 8 "Waste treatment".**
- Wash cells twice with 1 ml of pre-warmed PBS per well.
- Add 500 µl per well of complete cell culture medium containing the corresponding LPS and NM concentrations according to the pipetting scheme shown in Figure 3.
- Culture cells for appropriate time points (3 h, 8 h, 24 h) under standard growth conditions.

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		12/20




Figure 3 Application of stimuli. NMs as well as LPS are applied in 500 μl complete cell culture medium per well after two washing steps in PBS.

¹⁾ NM concentrations given here refer to metal oxide NMs. Carbon based NM concentrations are detailed in the text.

6.5 RNA isolation

Proceed according to the RNeasy® Micro Handbook provided by Qiagen. Only cell harvest is described briefly below.

- Prepare RLT-working buffer by adding 10 μl β -mercaptoethanol per 1 ml of RLT buffer. 350 μl RLT-working buffer is needed per well of a 24-well plate. For one plate as shown in figure 3 3.5 ml RLT buffer are mixed with 35 μl β -mercaptoethanol.
- After appropriate time points (3 h, 8 h, 24 h) remove supernatant (can be transferred to a separate 1.5 ml microreaction tube and used for protein measurement by ELISA. See SOP “Detection of cytokine release in THP-1 cells) and add 350 μl RLT-working buffer per well.
- Remove cells from the bottom of the cell culture well using a sterile pipet tip. The RLT-working buffer – cell mixture is viscous and clear.
- Make sure to scratch around the well long and rigorous enough to remove all cells from the well.
- Transfer everything into a 1.5 ml microreaction tube.
- Vortex to homogenize samples.
- Samples in RLT-working buffer can be frozen at -20°C till further processing.
- Proceed according to the RNeasy® Micro Handbook (attached to this SOP). For a brief overview see Figure 4.

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		13/20




Figure 4: Workflow of total RNA isolation. Picture source: RNeasy Micro Handbook 12/2007. (Modified to illustrate, that we start from cell material only, but that other biological sources (tissue, LMD samples) could also be processed.) To extract RNA stick closely to the protocol of the manufacturer.

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		14/20

6.6 RNA measurement

Here we describe only the computer assisted measurement using of the Nanodrop ND-1000 spectrophotometer. When using any other equipment please refer to the corresponding protocol/SOP.

Make sure to mix every sample properly before measurement. All samples are measured in duplicates.

- Start computer and nanodrop software.
- Clean sensor with ddH₂O and tissue.
- Choose “Nucleic Acid” as the measurement parameter.
- Initialization:
 - Put 1 µl ddH₂O on top of the sensor.
 - Close lid.
 - Press “o.k.” (software).
- Choose sample type “RNA40”.
- Clean sensor with ddH₂O and tissue.
- Blank measurement:
 - Put 1 µl of solvent (in our case RNA is dissolved in ddH₂O) on top of the sensor.
 - Close lid.
 - Press “blank” (software).
- Clean sensor with ddH₂O and tissue.
- Sample measurement:
 - Put 1 µl per sample on top of the sensor.
 - Close lid.
 - Press “measurement” (software).
 - Measure every sample twice.
 - Clean sensor with ddH₂O and tissue after each sample.
- Proceed with subsequent sample.
- Save report and transfer values to excel for further calculations.

Report will contain the following parameters:

- Concentration of RNA (ng/µl)
- A260: OD measured at 260 nm
(used to calculate RNA concentration according to equation (1))

$$RNA\ conc.\ (ng/\mu l) = OD(260\ nm) \times dilution\ factor \times 40 \quad (1)$$

- A280: OD measured at 280 nm
(used to assess contamination by proteins or phenols which absorb at 280 nm)
- Ratio 260/280
(used to assess the purity of the RNA preparation; a ratio ~2 is considered as “pure RNA”)

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		15/20

6.7 cDNA synthesis using iScript™ (BioRad)

Volumes given are for one single reaction only. Make sure to prepare enough iScript master mix for all samples to be processed.

iScript master mix:

- 5x iScript reaction mix 4 µl
- Add iScript reverse transcriptase 1 µl

Note: Prepare this mix shortly before usage. If necessary keep it on ice (at 4°C) until needed.

Note: As it is almost impossible to completely eliminate genomic DNA from RNA samples, it is important to run a so called “no amplification control” (NAC). This reaction contains all reagents except the reverse transcriptase. If later on in the qRT-PCR reaction a product is formed, it indicates the presence of contaminating DNA in the sample. Make sure to include these mock samples into the calculation for the iScript master mix.

Calculate the volumes needed per sample for a final amount of **200 ng RNA**.

All further steps are performed in a qRT-PCR 96-well plate.

- Add **200 ng RNA** of each sample into one 96-well.
- Ad up with **ddH₂O** to a final volume of **15 µl**.
- Add **5 µl iScript** master mix per well.
Results in a total reaction volume of **20 µl** (per well).
- Run reverse transcription program on iCycler:
5 minutes 25°C
30 minutes 42°C
5 min 85°C
hold at 4°C
- After reaction is completed add **80 µl ddH₂O** to the reaction mix.
Results in a total volume of 100 µl that contain an equivalent of 200 ng RNA. Assuming a complete transcription into cDNA this corresponds to a final concentration of **2 ng/µl cDNA**.

6.8 Real-time PCR using iQ™ SYBR® Green Supermix and iCycler

Volumes given are for one single reaction only. Working stock concentration for all primers used is 10 µM. All reactions are performed in duplicates in a qRT-PCR 96-well plate and in a total volume of 12 µl. An equivalent of 10 ng RNA (corresponds to 5 µl of the cDNA reaction performed in 6.7) per well is used.

Controls to be run:

- NAC (see 6.7): these wells receive 5 µl mock cDNA samples (as described above) instead of template. A signal here indicates contamination of samples with genomic DNA.

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		16/20

- No template control (NTC): these wells receive 5 µl ddH₂O instead of template. A signal here indicates that one or more of the qRT-PCR reagents is contaminated with the amplicon.

Actual procedure:

- Make sure to prepare enough **iQ SYBR Green master mix** for all samples to be processed:

Reagent	Volume
iQ SYBR Green Supermix	6.00 µl
primer_sense (10 µM)	0.24 µl
primer_antisense (10 µM)	0.24 µl
ddH ₂ O	0.52 µl
total volume	7.00 µl

- Add 5 µl cDNA template (2 ng/µl, resulting from cDNA synthesis reaction, see 6.7 above) per well of a 96-well qRT-PCR plate.

Note: This small volume has to be pipetted as a compact drop directly to the wall of the well.

- Add 7 µl of the iQ SYBR Green master mix. Avoid bubble formation.
- Seal the plate with an optically clear heat seal (e.g. BioRad). Avoid finger prints on top of the foil. These would impair fluorescence readout.
- Run appropriate PCR program. On the iCycler used here: *2Step60+Melt* (Figure 5)

step no.	temperature (°C)	duration (min.)
1	95	3:00
2	95	0:10
3	60	0:30
4	GO TO 2	39 more times
5	95	0:10
6 (melt curve)	65 to 95 (increment 0.5°C)	
END	4	∞




Figure 5: iCycler Program *2Step60+Melt*. 60 refers to the annealing temperature of the primers (60°C) in step 3.

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		17/20

6.9 Data evaluation

Resulting C_T values of the housekeeping gene (GAPDH) and the gene of interest (IL-8) are analyzed using the *ddCt* method. For a more detailed review on this algorithm please refer for example to Livak and Schmittgen, 2001 or Zhang et al., 2013.

ΔC_t refers to the difference between the C_T value of the gene of interest ($C_{T,x}$) and the C_T value of the housekeeping gene ($C_{T,0}$) under the same experimental conditions (see equation 2).

$$\Delta C_{t_x} = C_{T,x} - C_{T,0} \quad (2)$$

$\Delta\Delta C_t$ represents the difference between the ΔC_{t_x} values of one treatment condition (e.g. TNF- α treatment to induce gene expression) and the ΔC_{t_y} values of a reference treatment condition (e.g. untreated samples) (see equation 3).

$$\Delta\Delta C_{t_x} = \Delta C_{t_x} - \Delta C_{t_y} \quad (3)$$

The relative expression of the gene of interest in treated vs. reference samples is finally given by equation (4):

$$relative\ expression(x) = 2^{-\Delta\Delta C_{t_x}} \quad (4)$$

Note: For the *ddCt* calculation to be valid, the amplification efficiencies of the gene of interest and the housekeeping gene must be approximately equal. A protocol how to assess the efficiency of a specific primer pair can be found for example at

http://www3.appliedbiosystems.com/cms/groups/mcb_marketing/documents/generaldocuments/cms_053906.pdf or in brief in the related document “qRT-PCR_efficiency”.

7 Quality Control, Quality Assurance, Acceptance Criteria

RNA purity: A ratio of absorbance at 260 nm and 280 nm of ~ 2 is generally accepted.

PCR efficiency: The PCR efficiency of each primer pair has to be determined in advance and is not part of this SOP. Generally, a PCR efficiency between 90 and 110% is considered acceptable (*Lifetechnologie*).

Negative controls (NAC and NTC) have to be negative (no C_T value detectable).

Melt curves of each amplicon yield a sharp peak at the melting temperature of the respective amplicon. This indicates that the products are specific. As a conclusion SYBR Green fluorescence is directly correlated to the accumulation of the product of interest.

Biological readout: Treatment of THP-1 cells with 1 $\mu\text{g/ml}$ LPS for 3 hours should result in an approximately 50-fold induction of IL-8 expression, 50-fold induction of TNF- α expression and approximately 3-fold induction of Nrf2 expression compared to the untreated 3 hour control sample.

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		18/20

8 Health and Safety Warnings, Cautions and Waste Treatment

Cell seeding has to be carried out under sterile conditions in a laminar flow cabinet (biological hazard standard). For this only sterile equipment must be used and operators should wear laboratory coat and gloves (according to laboratory internal standards). Special care has to be taken during PMA handling (carcinogenic potential of the substance!).

PMA waste treatment: use a separate exhaust extraction system with a collecting flask containing already 20 ml 5 M NaOH to neutralize PMA. The resulting non-toxic solution is very alkaline and has to be neutralized using HCl before final disposal in the sink.

NaOH is **corrosive**. It causes severe burns. Wear especially eye/face protection. Dissolution of NaOH is an exothermic reaction, the solution will get fairly hot – be careful! It is strongly recommended to wear eye protection when handling 5 M NaOH.

HCl is **corrosive and irritant**. It is very hazardous in case of skin contact, of eye contact and of ingestion. It is slightly hazardous in case of inhalation. Therefore avoid inhalation as well as contact with skin and eyes and avoid exposure in general.

Discard all materials used to handle cells (including remaining cells themselves) according to the appropriate procedure for special biological waste (i.e. by autoclaving).

9 Abbreviations

cDNA	copy deoxyribonucleic acid
ddH ₂ O	double-distilled water
DMSO	dimethyl sulfoxide
DNA	deoxyribonucleic acid
ELISA	enzyme-linked immunosorbent assay
FCS	fetal calf serum
g	constant of gravitation
GAPDH	glyceraldehyde 3-phosphate dehydrogenase
HO-1	heme oxygenase 1
IL-8	interleukin 8
LPS	lipopolysaccharide
mRNA	messenger ribonucleic acid
NAC	no amplification control
NM	nanomaterial
Nrf2	nuclear factor erythroid 2 (NF-E2)-related factor 2
NTC	no template control
ON	overnight
PBS	phosphate buffered saline
PMA	phorbol 12-myristate 13-acetate
ppm	Parts per million
PSN	Penicillin, Streptomycin, Neomycin
qRT-PCR	quantitative real-time reverse transcription polymerase chain reaction
ROS	reactive oxygen species
RPMI	Roswell Park Memorial Institute medium
TNF- α	tumor necrosis factor alpha

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		19/20

10 References

Donaldson K, Stone V, Tran CL, Kreyling W, Borm PJ; 2004: Nanotoxicology. *Occup Environ Med* 61(9):727-728

Knippers R.: Molekulare Genetik. Textbook. Thieme Verlag.

Limbach LK, Wick P, Manser P, Grass RN, Bruinink A, Stark WJ; 2007: Exposure of engineered nanoparticles to human lung epithelial cells: influence of chemical composition and catalytic activity on oxidative stress. *Environ Sci Technol* 41 (11):4158-4163

Livak KJ, Schmittgen TD; 2001: Analysis of Realtime Gene Expression Data Using Real-Time Quantitative PCR and the $2^{-\Delta\Delta CT}$ Method. *Method*, 25(4): 402-408

Schins RP and Knaapen AM; 2007: Genotoxicity of poorly soluble particles. *Inhal Toxicol* 19 (Suppl 1):189-198

Zhang JD, Ruschhaupt M, Biczok R;2013: ddCt method for qRT-PCR data analysis.
<http://www.bioconductor.org/packages/release/bioc/vignettes/ddCt/inst/doc/rtPCR.pdf>

Document Type	Document ID	Version	Status	Page
SOP	I+O_qRT-PCR_THP-1	1.0		20/20